


Disposable Hypodermic EMG Needle Electrode



Technomed understands the importance of a crystal clear signal, ease of use, repositionability and patient comfort.

Technomed's Hypodermic Needle Electrode is a high quality product developed for EMG guided drug delivery.

Technomed's Hypodermic Needle Electrode can be used for multiple insertions without degradation of its characteristics.


The hub is designed to give a good and secure feel in the hand.

User Benefits

- Optimal patient comfort with ultra sharp tip and special coating for very smooth penetration and repositioning
- Leak-free and minimum drug residual volume
- Optimized to fit to all slip-tip and Luer Lock syringe connections
- Efficiency of use with color coding
- Easy to open pouch for quick and hassle free preparation of the procedure

Available variants

Needle Length	Diameter	Lead Wire	Hub color	Packing	Part Number
25 mm	0.30 mm (30 Gauge)	0.75 m	Red	10 individually packed electrodes per box	TE/L2530-335
25 mm	0.40 mm (27 Gauge)	0.75 m	Yellow	10 individually packed electrodes per box	TE/L2540-335
37 mm	0.40 mm (27 Gauge)	0.75 m	Orange	10 individually packed electrodes per box	TE/L3740-335
37 mm	0.45 mm (26 Gauge)	0.75 m	Green	10 individually packed electrodes per box	TE/L3745-335
50 mm	0.50 mm (25 Gauge)	0.75 m	Blue	10 individually packed electrodes per box	TE/L5050-335
75 mm	0.70 mm (22 Gauge)	0.75 m	Violet	10 individually packed electrodes per box	TE/L7570-335



Specifications

Lead-free	✓
Latex-free	✓
BPA-free	✓
DEHP < 0.1% by weight	✓
REACH & RoHS compliant	✓
EtO sterilized	✓
Lead wire core	Tin Plated Copper
Lead wire insulation	PVC
Connector	DIN 42802
Shelf life for unopened pouches	36 months

Your Technomed Business Partner:

USA
Rx ONLY

CE
0344



Technomed
medical accessories

Taking you further. Step by step.

Technomed

Amerikalaan 71
6199 AE | Maastricht-Airport | The Netherlands
T +31 43 408 68 68 | info@technomed.nl

Technomed USA

Suite 4990 Constellation Drive
White Bear Township | MN 55127 | United States
T +1 651 288 0890 | info@technomed-usa.com

www.technomed.nl